



1. Cellular Indicator shows cellular signal strength.
2. Battery Indicator, which surrounds the Call Button, shows level of battery charge.
3. Call Button, when pushed, initiates a two-way call to the response team.
4. Speaker allows the end user to hear the response team.
5. Microphone allows the response team to hear the end user.

Cellular Indicator

The Cellular Indicator on top of Libris shows cellular signal strength.

- When the Cellular Indicator blinks **GREEN**, Libris has a strong signal.
- When the Cellular Indicator blinks **AMBER**, Libris has a moderate signal.
- When the Cellular Indicator blinks **RED**, Libris has no signal. Move to an area with cellular coverage.

Battery Status

The Libris Battery Indicator, which surrounds the Call Button, shows the level of battery charge.

- When the Battery Indicator flashes **GREEN**, Libris is well charged.
- When the Battery Indicator flashes **AMBER**, charge Libris soon.
- When the Battery Indicator flashes **RED**, the battery is low. Charge Libris immediately.

Libris also plays two audio messages when battery is running low:

1. "Charge your personal mobile responder soon."
2. "Charge your personal mobile responder NOW. Please place your personal mobile responder in the charger now."